



# Updating ImageCentrik

## Method 1


1. Open your ImageCentrik program and navigate to Help > Update



2. Click Next to continue.




3. Click Next to continue




4. Close ImageCentrik by clicking the red X box in the upper-right hand corner of the screen

5. Click Yes to continue



Setup will now install the Update for ImageCentrik!

6. If you run ImageCentrik Network version a window will pop up to display any other requirements needed (like updated server software). Click Ok and then finish.



## Method 2

1. Contact [support@dlcenters.com](mailto:support@dlcenters.com) to receive a link to the latest version of ImageCentrik

*Note: Do not trust ImageCentrik downloads that are not from [www.DLCenters.com](http://www.DLCenters.com) or its staff*

2. Double-click the download executable.
3. Continue through the setup process as if you were to install the program from scratch

